


Vad är UriCap Female?

UriCap är en innovativ och säker urinuppsamlingslösning för kvinnor. Ett alternativ till urinkateter och inkontinensskydd.



Enkel läckagefri engångsanvändning i upp till 24 timmar

Tillverkad av mjuk medicinsk silikon

Fixering genom naturlig adhesion och gravitation, kräver inte lim eller tejp

Försluter tätt runt uretramynningen och ansluts till en urinuppsamlingspåse

Varför UriCap Female?

- Säkert alternativ till kateterisering
- Snabb och enkel applicering och borttagning
- Tillhandahåller ett rent urinprov
- Kan minska antalet personalkrävande byten av inkontinensskydd
- Ingen steril procedur behövs
- Reducerar avfallet och minimerar lagringsutrymmet
- Alternativ vid IAD (Incontinence Associated Dermatitis)

Product datasheet Uricap Female




MDD Class I

Components:

The UriCap Female is assembled from two parts, the UriCap Female and the Extension Tube.

UriCap Female Device Specifications

Extension Tube SX Specifications



Name	UriCap Female Device	Name	Extension Tube SX-L
PN	FPR02020L	PN	FPR00004
Material	Medical grade silicone	Material	Tube - Medical grade silicone Connector- Medical grade polypropylene
Size packed	15*10*2.5 cm.	Size packed	25*7.5*1.5 cm.
Weight packed	7 gr.	Tube Length	28 cm.
Shelf life	3 years	Weight packed	9.6 gr.
Storage conditions	10-30°C	Shelf life	3 years
Comments	Single Use device. certified for maximum 24 hrs. continuous use.	Storage conditions	10-30°C
		Comments	.

Contact in Sweden

Invent Medic Sweden
Nytänkargatan 4
223 63 Lund
info@inventmedic.com

Manufacturer

TillaCare Ltd.
2 Etgar st., 390321 Tirat
Carmel, Israel
E-mail: sales@tillacare.com

